

Hill-Rom Advanta™ Bed Systems



At a glance

- Retractable frame
- 4 wheel dual-locking brakes
- Low height of 15.25"
- OneStep® instant CPR activation
- One-button cardiac chair
- Ergonomic head and footboard
- Line-of-Sight® head angle indicator
- Addresses six main entrapment areas

Engineered for safety. Designed for comfort.

Hill-Rom's Advanta bed is engineered to help reduce the likelihood of patient falls and offers the added benefit of an optional four-zone pressure reduction surface. While the Advanta bed was created to help address your patients' needs, it also offers features that contribute to the safety and efficiency of caregivers.

Hill-Rom Advanta™ Bed Systems



Features

- **Patient Position Monitor (PPM)** - Alarms may be customized to alert staff when the patient is moving, exiting or out of bed. Because the nurse selects the sensitivity of the alarms, there are fewer “false” calls on anticipated patient movements. Patients at risk can quickly be identified and assessed.
- **Low bed height** - A 15.25” bed height allows patients to place both feet firmly on the ground as they move in and out of bed. Based on their capabilities, patients may require less assistance from staff.
- **Dual-locking Brake System** - Designed to help prevent the bed from moving during egress and ingress, the braking system also allows the bed to aid in supporting the patient's weight should he or she become unsteady when moving in or out of bed.



Options

- **Optional ZoneAire® sleep surface** offers pressure relief for at-risk patients.
- **Point-of-Care® siderails** are conveniently located for caregiver efficiency.
- **In-bed scale controls** optional scale weighs patients without having to move them. The Advanta bed's scale can weigh in any bed position except when the bed is in full Trendelenburg while also in the lowest position.

Technical data

Sleep surface	36" x 84" (91.6 cm x 213.4 cm)
Frame length	(flat): 91" (231.1 cm), (retracted): 78.5" (199.4 cm)
Frame width.....(rails up):	42.75" (108.6 cm), (rails down): 40.25" (102.2 cm)
Underbed clearance (with 5" casters)	4.75" (12.1 cm)
High position (floor to sleep deck)	28.5" (72.4 cm)
Low position (floor to sleep deck)	15.25" (38.7 cm)
Maximum head elevation65°
Maximum Trend/Reverse Trend	+/- 11°
Caster size	5" (12.7 cm)
Maximum patient weight limit	500 lbs. (226.76 kg)
Scale accuracy:	1% of patient weight
Scale repeatability	0.2% of patient weight
Scale weight capacity	450 lbs. (181.4 kg)
Electrical leakage	<65 microamps @ 115VAC



Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2003 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.

CTG471rb 10/06/03