



Kangaroo™
Enteral Feeding Pumps



Continuum of Care

The Kangaroo™ Family of Enteral Feeding Pumps

For the past thirty years, Covidien has offered market leading and innovative Kangaroo feeding pumps. The current Kangaroo platforms continue the tradition of excellence, offering a full spectrum of enteral feeding capabilities including pumps designed to feed and flush, an ambulatory pump for patients on the go, and a standard, dependable, easy-to-use pump.

In addition to providing innovative solutions for all patient populations, Covidien offers a complete nutritional delivery system that provides both engineered controls and a visual indicator to help prevent tubing misconnections.

For more information on our nutritional delivery platform visit www.SafeEnteralConnections.com

Kangaroo™ 924

Kangaroo ePump™

Kangaroo Joey™



Kangaroo ePump™ Enteral Feeding Pump

The right pump for acute and long term care centers

Revolutionizing the delivery of enteral feeding begins with e

When patients require nutrition through continuous feeding, intermittent feeding or feeding and flushing, the Kangaroo ePump feeding pump delivers in one compact, easy-to-use device.

Kangaroo Joey™ Enteral Feeding Pump

The right pump at home

Change the world of ambulatory feeding and improve patient outcomes with the pump that feeds and hydrates

An innovation from Covidien that provides feeding and programmable flushing in a completely portable device. The first hydrating ambulatory enteral device available, Covidien's Kangaroo Joey pump assures clinicians of accurate and reliable function while providing patients with freedom and mobility.

Kangaroo™ 924 Enteral Feeding Pump

The right pump made easy

The standard in Enteral Feeding

When ease of use, accuracy and dependability matter, the Kangaroo 924 Enteral Feeding Pump delivers it all.





Kangaroo ePump™ Enteral Feeding Pump

The right pump for acute and long term care centers

ATM-style Interface

Makes programming of the pump, intuitive, easy and quick, saving nursing time and enhancing proper programming inputs.

Programmable Flush

A one-of-a-kind option that allows programming specific flushing rates (10-500 mL) and flushing intervals (1-24 hrs), saving nursing time and eliminating the need for most manual flushes.

Attitude Independent

Allows for patient portability and ensures uninterrupted feeding and hydration should the pump accidentally be placed on its side, front or back.

Auto Prime

Reduces the need for time consuming manual priming

72 Hour Feeding History

Stores the previous 72 hours of feeding and flushing history, ensuring the patient has received the prescribed dose.





Kangaroo Joey™ Enteral Feeding Pump

The right pump at home

Shared Features with the Kangaroo ePump Enteral Feeding Pump

Enable the patient on the go with an ambulatory enteral device which has market-leading, innovative features, including an ATM-style interface, auto priming, attitude independence, and a 72-hour feeding history, in addition to:

Programmable Flush

First ambulatory pump on the market with a programmable flush feature that allows programming specific flushing rates (10-500 mL) and flushing intervals (1-24 hrs), delivering nutrition and hydration while patients are on the go.

Enhanced Accuracy and Extended Battery Life

+/- 7% with an 18 hour battery life.

Shock and Water Resistant

Pump is constructed with a durable rubber molding improving durability and protecting against ingress fluids (IPX4 protection).

Fully Portable Design with Backpack and Quick-Release Pole Clamp

Allows for convenient, discrete use.



Kangaroo™ 924 Enteral Feeding Pump

The right pump made easy

Intuitive, User Friendly Controls

Front panel button controls are easy to read and intuitive to program and operate.

Over-infusion Safeguard

The Volume To Be Delivered (VTBD) feature ensures correct and consistent delivery of nutrition while enhancing patient safety

Unique Safety Sensor

Optical sensor provides safeguards against improper placement of the pump set, which results in a higher level of patient safety and product performance.

24 Hour Memory

Pump memory will retain settings and volume delivered even if pump is turned off, helping to facilitate keeping records of formula infused.

Audible and Visual Alarms

Alarms will reinforce clinician or caregiver confidence in the safe delivery of enteral nutrition.

Kangaroo ePump™ Enteral Feeding Pump

Code Number	Description	Case Quantity
382400	Kangaroo ePump with Pole Clamp	1
382491	Power Adapter	1
382492	Pole Clamp	1
F010484WT	Battery	1

PUMP SETS

774655	Proximal Spike Set	30
772055	500 mL Pump Set	30
773656	1000 mL Pump Set	30
776255	600 mL Rigid Container Pump Set	30
776250	1200 mL Rigid Container Pump Set	30
774669	Proximal Spike Feed and 1000 mL Flush Set	30
773662	1000 mL Feed Set w/ 1000 mL Flush Set	30
672062	Spike Feed and 500 mL Flush Set	30
674668	500 mL Feed Set with 500 mL Flush Set	30
714650	Sterile Proximal Spike Set	30
712054	Sterile 500 mL Pump Set	30
713650	Sterile 1000 mL Pump Set	30
716154	Sterile 100 mL Burette Set	30

Kangaroo™ 924 Enteral Feeding Pump

Code Number	Description	Case Quantity
392455	Kangaroo 924 Enteral Feeding Pump with Pole Clamp	1
773621	1000 mL Bag Set (Indiv. Pack)	30
703621	1000 mL Bag Set (Bulk Pack)	30
774625	Spike Set (Indiv. Pack)	30
772025	500 mL Bag Set	30



Kangaroo Joey™ Enteral Feeding Pump

Code Number	Description	Case Quantity
383400	Kangaroo Joey Enteral Feeding Pump w/ Pole Clamp	1
383493	Kangaroo Joey Pole Clamp	1
383494	Kangaroo Joey Flexible Pole Clamp	1
383491	Kangaroo Joey Power Cord	1
F010506	Kangaroo Joey Battery	1
770026	Blue Mini Backpack	1
770028	Blue Super Mini Backpack	1
770034	Pink Mini Backpack	1
770029	Pink Super Mini Backpack	1
770027	Green Mini Backpack	1
770032	Green Super Mini Backpack	1
770025	Black Mini Backpack	1
770031	Black Super Mini Backpack	1
PUMP SETS		
764655	Kangaroo Joey Proximal Spike Set	30
762055	Kangaroo Joey 500 mL Pump Set	30
763656	Kangaroo Joey 1000 mL Pump Set	30
764669	Kangaroo Joey Spike Set w/ 1000 mL Flush Bag	30
763662	Kangaroo Joey 1000 mL Feed Set with 1000 mL Flush Set	30

Covidien is pleased to provide a complete nutritional delivery system that provides both engineered controls and a visual indicator to help prevent tubing misconnections. The following accessories are designed to be compatible with the Nestlé®* SpikeRight™* system.

Kangaroo™ Feeding Sets for use with Nestlé®* SpikeRight™* System

Code Number	Description	Case Quantity
775655	Kangaroo Safety Screw Spike for Kangaroo ePump	30
765555	Kangaroo Joey Safety Screw Spike	30
775662	Kangaroo Safety Screw Spike w/1000 mL Flush Bag for Kangaroo ePump	30
765662	Kangaroo Joey Safety Screw Spike w/1000 mL Flush Bag	30
775756	Kangaroo Safety Screw Spike for Kangaroo 924 Pump	30
775762	Safety Screw Spike Adapter Cap	30



For more information, please visit
www.SafeEnteralConnections.com
or call 1 800 962-9888



COVIDIEN

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other TM marked brands are trademarks of a Covidien company.

®*TM* Nestlé® and SpikeRightTM are trademarks of the Société des Produits Nestlé S.A., Vevey, Switzerland.

© 2009 Covidien.

H6450-XM-0809

15 HAMPSHIRE STREET
MANSFIELD, MA
02048

1-800-962-9888
508-261-8000

WWW.COVIDIEN.COM